



LITHOTOMY STIRRUP

- * Next-Level Innovation for Next-Level Care
- * Dual-Bar Design: Help prevent medial leg drop with this innovative safety feature
- * Boot Release Handle: Easily adjust leg position with a single-point latch
- * Welded Pads, Removable Silicone Straps: Control potential infection with easy-to-clean welded pads and silicone straps
- * Patient positioning using the Lithotomy Stirrups features an extended lateral boot design that supports the patient's leg entirely, preventing superficial nerve damage and prevents unwanted leg rotation intra-operatively.

TECHNICAL PARAMETER

ITEM	PARAMETER
Cantilever Length	1000mm
Back section up/down	90° / -35°
Lift up and down	30° / 15°
Telescopic distance of sheath	330mm
Patient Weight Capacity	150Kgs
Specification of guardrail	25*10mm /25*8mm /30*10mm /30*8mm /28.5*9.5mm
Product Weight	13kgs



SUPPORT LEVER CONTROL HANDLE



FIXED HANDLE



GERMAN PNEUMATIC AUXILIARY DEVICE



CLAMP



REUSABLE PAD



YA-A601

No.3 Jinxiu Road, Jinfeng Industrial Zone, Jinfeng Town, Zhangjiagang, Jiangsu Province, China

+86 13962231292
info@medikmedical.com

